

MediCap™ USB170

Medical Still Image Recorder




Record images the fast and easy way

The MediCap USB170 is the easy way to capture still images of your medical procedures. Images are recorded directly to a USB flash drive or optional hard drive. When you finish recording, just plug the drive into any computer – then you can view, archive, edit, share, etc. It's so simple to use, you'll be saving images just minutes after unpacking your USB170. Now it's easy to get started with digital imaging, and at price that's remarkably affordable.

Versatile and affordable...

- ▶ Works with endoscopes, arthroscopes, ultrasounds, or any video source with standard outputs
- ▶ Store thousands of images on the included USB flash drive or optional USB hard drive
- ▶ High-quality images in JPG, TIF, PNG, or DICOM
- ▶ Collect patient information using any USB keyboard
- ▶ Record via button, footswitch, or trigger
- ▶ Network connection option
- ▶ Compact, rugged design fits easily on any cart
- ▶ No training required – a snap to connect and use



MediCap™ USB170 Specifications

Device Type	Medical Still Image Recorder Records digital still images from analog video sources to USB flash drive or optional USB hard drive
Case	Rugged metal case with high-impact plastic faceplate
LCD Monitor	2 1/2", full color (for navigation only, not for diagnosis)
Buttons	Sealed membrane, fluid resistant
Video In/Out	PAL or NTSC (menu selectable) S-video or composite BNC connections
USB Support	USB2.0, high speed, 1 front connector, 3 rear connectors (all powered)
Networking	Ethernet T10/100 connection via Network Upgrade Kit (optional)
Footswitch	Hands-free capture Standard 3.5mm mini jack
Mechanical	Size: 240 x 200 x 63mm / 9.5 x 7.9 x 2.5" Weight: 1.7kg / 3lb 12oz (appx.)
Environmental	Storage: -40° to +85° C Operating: -20° to +40° C
Power	Power 110/220V 3-prong AC jack
Media	USB flash drives or external USB hard drive
Image Formats	JPEG, TIFF, PNG, DICOM
Image Size	1024 x 768, 800 x 600, 640 x 480 pixels
Certifications	UL Medical, CE, FCC, FDA Class 1, RoHS
Kit Contents	MediCap USB200 USB flash drive Video cables Power cable User guide

Also Available...



MediCap USB200

All the features of the USB170, *plus* video recording. Holds two hours of video clips on the included flash drive or hundreds of hours on an optional hard drive.